


TS-37C(NIX-521)

適用遠心機	NIX-521
バケット架け数	4箇所／ローター

バケット名	ラック容量 / ラック名	最高 回転数 [rpm]	最大遠心 加速度 [G]	回転半径 [mm]		備考
				Rmax	Rmin	
B437 	250 ml チューブ × 1本 3625C-01P 	3,200	2,070	181	106	
	100 ml チューブ × 2本 3610C-G02P 	3,200	2,050	179	106	
	50 ml チューブ × 5本 3650-G05P 	3,200	2,030	177	106	
	50 ml エカチューブ × 5本 15 ml エカチューブ × 2本 3650-TC05P 	3,200	2,100	183	106	※1
	15 ml チューブ × 16本 3615-G16P 	3,200	2,080	182	106	
	15 ml カチャチューブ × 14本 3615-TC14P 	3,200	2,100	183	106	
	14ml チューブ × 18本 3614C-18P 	3,200	2,020	176	118	
	6 ml チューブ × 35本 3606C-35P 	3,200	2,100	183	124	
	5 ml チューブ × 48本 3605C-48P 	3,200	2,100	183	130	
	2 ml マイクロチューブ × 36本 3602C-36P 	3,200	1,920	168	130	
	ティーフウエルプレート × 2枚 ※2 マイクロプレート × 4枚 ※3 AS36C-96D 	3,200	2,080	182	—	※2 ※3
	50 ml チューブ × 6本 3750-G06P 	3,200	2,030	177	106	
	50 ml エカチューブ × 6本 15 ml エカチューブ × 2本 3750-TC06P 	3,200	2,100	183	106	※4
	15ml チューブ × 24本 3715-G24P 	3,200	2,080	182	106	
	15 ml エカチューブ × 16本 3715-TC16P 	3,200	2,100	183	106	

※1 50 ml × 4本と15 ml × 2本までなら混載可能。 ※2 負荷寸法 (L×D×H) 128.5×86×90mmまで。

※3 負荷寸法 (L×D×H) 128.5×86×58mmまで。 ※4 50 mlと15 mlチューブは混載不可。

TS-37C 使用チューブ一覧 (NIX-521)

バケット	ラック	アダプター	チューブ						BART コード	実 容量 [ml]	最大本数 (本×箇所)	許容 回転数 [rpm]	最大遠心 加速度 [G] (半径 [mm])	備 考		
			公称 容量 [ml]	メーカ 略称	型名・名称 (材質)	底 形	寸法 [φ×Lmm]	許容遠心 加速度 [G]								
B437 	3625C-01P 	3124-0010 		NAL	3123-0250 ※1 (PC)	R	φ 61.9×145.8	27,500		250			2,040 (178)	※1		
				NAL	3120-0250 (PPCO) 3121-0250 (HDPE) 3122-0250 (PC) 3127-0250 ※2 (FEP)	F F F F	φ 61.8×127.7 φ 61.2×131.9 φ 61.8×127.6 φ 60×128.8	13,200 8,000 27,500 4,000					※2			
		352090 		FLC	352075 (PP)	C	φ 61×137	7,500	1	225	1×4	3,200	2,070 (181)	※3		
				FLC	352076 (PP)	C	φ 61×118	7,500		175						
		3610C-G02P 		NEG	丸形遠心管 (G)	R	φ 45×120	—	2	100	2×4	3,200	2,050 (179)			
		3614C-18P 		FLC	352057 (PS)	R	φ 17×100	1,400	7	14	18×4	2,600	1,330 (176)			
					352059 (PP)	R	φ 17×100	3,000				3,200	2,020(176)			
					352017 (PS)	R	φ 17×100	1,400				2,600	1,330 (176)			
					352018 (PP)	R	φ 17×100	3,000				3,200	2,020(176)			
		3606C-35P 		EIK	JA1000 (PS) JA2000 (PS) JA3000 (PS) JA4000 (PS)	R	φ 14.5×83	—	8	6	35×4	3,200	2,100 (183)			
	3605C-48P 		FLC	352002 (PP)	R	φ 12×75	3,000	9	5	48×4	3,200	2,100 (183)	※4			
				352053 (PP)	R	φ 12×75	3,000			24×4						
				352063 ※4 (PP)	R	φ 12×75	3,000									
				352058 ※4 (PS)	R	φ 12×75	1,400									
				352235 ※4 (PS)	R	φ 12×75	1,400							48×4	2,600	1,380 (183)
				352008 (PS)	R	φ 12×75	1,400									
				352052 (PS)	R	φ 12×75	1,400									
	3602C-36P 		SIO	SHIONOGI チューブ (PS)	R	φ 12×72	—		4	48×4	3,200	2,100 (183)				
			EPP	セイワロックチューブ (PP)	C	φ 11×42	25,000	A	2	36×4	3,200	1,920 (168)				
			EPP	セイワロックチューブ (PP)	C	φ 11×41.5	25,000		1.5							
	マイクロチューブ 3810 (PP)			C	φ 11×41.5	25,000										
	AS36C-96D 	—	—	※5 ディープウェルプレート	—	—	—	B		2枚×4	—	— (—)	※5			
				※5 マイクロプレート						4枚×4						

※1 アダプター (3124-0010) が別途必要。 ※2 4℃設定のみ使用可。 ※3 コニカルチューブ用アダプター (352090) が別途必要。 ※4 ラックの穴を1つおきに使用してください。

※5 使用するプレートが遠心に耐えられるかどうか事前にご確認ください。

● イラストは商品イメージであり、実際の色合い、形状など異なる場合があります。

TS-37C 使用チューブ一覧 (NIX-521)

バケット	ラック	アダプター	チューブ						BART コード	実 容量 [ml]	最大本数 (本×箇所)	許容 回転数 [rpm]	最大遠心 加速度 [G] (半径 [mm])	備 考	
			公称 容量 [ml]	メーカ 略称	型名・名称 (材質)	底 形	寸法 [φ×Lmm]	許容遠心 加速度 [G]							
B437 	3650-TC05P 	 50ml	COR	430291	(PP)	C	φ 29×116	15,500	4	50	5×4	3,200	2,100 (183)		
				430304	(PET)	C	φ 29×116	3,600							
				430829	(PP)	C	φ 29×116	15,500							
				352070	(PP)	C	φ 30×115	9,400							
		 15ml	COR	430053	(PET)	C	φ 16×120	3,600	4	15	2×4	3,200	2,100 (183)		
				430766	(PP)	C	φ 16×120	12,000							
				430791	(PP)	C	φ 16×120	12,000							
				352196	(PP)	C	φ 17×120	6,000							
	FLC	352095	(PS)	C	φ 17×120	1,800		2,900	1,720 (183)						
		3650-G05P 	 50ml	NEG	丸形遠沈管 (G)		R	φ 35×100	—	3	50	5×4	3,200	2,030 (177)	
	3615-G16P 	 15ml	NEG	P-16.5S (G)		R	φ 16.5×105	—	5	15	16×4	3,200	2,080 (182)		
	3615-TC14P 	 15ml	COR	430053	(PET)	C	φ 16×120	3,600	6	15	14×4	3,200	2,100 (183)		
				430766	(PP)	C	φ 16×120	12,000							
430791				(PP)	C	φ 16×120	12,000								
FLC			352196	(PP)	C	φ 17×120	6,000	2,900				1,720 (183)			
			352095	(PS)	C	φ 17×120	1,800								
B437 	3750-TC06P 	 50ml	COR	430291	(PP)	C	φ 29×116	15,500	4	50	6×4	3,200	2,100 (183)		
				430304	(PET)	C	φ 29×116	3,600							
				430829	(PP)	C	φ 29×116	15,500							
				352070	(PP)	C	φ 30×115	9,400							
		 15ml	COR	430053	(PET)	C	φ 16×120	3,600	4	15	2×4	3,200	2,100 (183)		
				430766	(PP)	C	φ 16×120	12,000							
				430791	(PP)	C	φ 16×120	12,000							
				352196	(PP)	C	φ 17×120	6,000							
	FLC	352095	(PS)	C	φ 17×120	1,800		2,900	1,720 (183)						
		3750-G06P 	 50ml	NEG	丸形遠沈管 (G)		R	φ 35×100	—	3	50	6×4	3,200	2,030 (177)	
	3715-G24P 	 15ml	NEG	P-16.5S (G)		R	φ 16.5×105	—	5	15	24×4	3,200	2,080 (182)		
	3715-TC16P 	 15ml	COR	430053	(PET)	C	φ 16×120	3,600	6	15	16×4	3,200	2,100 (183)		
				430766	(PP)	C	φ 16×120	12,000							
430791				(PP)	C	φ 16×120	12,000								
FLC			352196	(PP)	C	φ 17×120	6,000	2,900				1,720 (183)			
			352095	(PS)	C	φ 17×120	1,800								

● イラストは商品イメージであり、実際の色合い、形状など異なる場合があります。