

遠心機

型式名称	消費電力 (W)	発熱量(kcal/h)
AX-310/311	930	800
AX-320/321	1,030	890
AX-501	760	654
AX-511	960	826
AX-521	970	834
CAX-370/371	1,040	900
CAX-571	990	851
Cryo Spin	8	7
EIX-136	1,070	920
EX-126	950	820
EX-136	1,030	890
GRX-220	3,860	3,320
GRX-250	3,860	3,320
KINTARO-18	210	180
KINTARO-24	300	260
KITMAN-18	350	300
KITMAN-24	380	330
KITMAN-T24	400	344
LC-100	600	520
LC-120	600	520
LC-121	600	520
LC-122	600	520
LC-200	470	400

製品消費電力／発熱量 一覧



LC-220	600	520
LC-230	730	630
LCX-100	430	370
LCX-200	660	568
LIX-130	1,070	920
LIX-135	1,150	990
LT-015/LT-010M	30	26
LX-120	860	740
LX-121	950	820
LX-130	1,030	890
LX-131	1,030	890
LX-140	1,030	890
LX-141	1,150	990
MDX-110	760	654
MDX-210	810	697
MDX-310	810	697
Multi Spin	8	7
MX-100	465	400
MX-105/107	650	560
MX-150	480	415
MX-160	600	520
MX-200	740	640
MX-201	790	680
MX-205/207	780	670
MX-300/301	860	740

製品消費電力／発熱量 一覧



MX-305/307	920	790
NIX-521	1,050	903
One Spin	4	3
PMC-060	8	7
PS-020	40	35
SRX-201	1,950	1,677
SUPREMA21	2,600	2,240
SUPREMA23	3,300	2,840
SUPREMA25	3,300	2,840
TX-160	600	520
TX-201	750	645

オートクレーブ

型式名称	消費電力 (W)	発熱量(kcal/h)
BS-235/245	1,500	1,290
BS-305/325	2,000	1,720
ES-215	1,500	1,290
ES-315	2,000	1,720
FLS-1000	3,500	3,010
KS-243	1,500	1,290
KS-323	2,000	1,720
LBS-245	1,500	1,290
LBS-325	2,000	1,720
LSX/SX-300	1,500	1,290
LSX/SX/BSX-500	2,000	1,720
LSX/SX-700	3,000	2,580

製品消費電力／発熱量 一覧



MSS-325	2,000	1,720
SD/SDL-320	2,000	1,720
SD/SDL-321	2,000	1,720
SR-240	3,000	2,580

その他

型式名称	消費電力 (W)	発熱量(kcal/h)
CC-105	500	430
CR-300L	300	160
CR-400L	400	210
CR-620/420	320	280
MS-100	460	400
MS-100R	600	520
MT-360/MT-400	10	8.6
MV-100	290	250
TCC-1207/1607(50/60Hz)	730	1,130/1,300
TCC-2407(50/60Hz)	830	1,440/1,590
TU-450	400	344
TU-500	300	258
TU-1000	600	516
UD-100	90	77
UD-201	300	258
UD-211	330	280
UDS-200	200	172
UR-20P/21P	36	30

人工気象器

型式名称	消費電力 (W)	発熱量(kcal/h)
CF/CFH-305	2,300	1,978
CF/CFH-405	2,700	2,320
CF/CFH-415	2,700	2,320
CL/CLH-301	2,300	1,978
CLE-303	1,250	1,075
CLE-305	1,150	989
CLE-405	1,250	1,075